Tierwaage
Veterinary scale
Balance pour animaux
Balanzas de animales
Scaglie di animali

6958

Gebrauchsanweisung	Seite 2
Operating Instructions	Page 24
Mode d'emploi	Page 46
Instrucciones de manejo	Página 72
Istruzioni d'uso	Pagina 94

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Thank you for choosing this Soehnle Professional product.
Soehnle Professional is a brand of Soehnle Industrial Solutions GmbH.
This product is equipped with all features of the latest technology.

If you have any questions or if you have problems with your device that are not covered in the operating instructions, please contact your Soehnle Industrial Solutions service centre or our customer service.

Telephone: +49 7191 3453-220
Fax: +49 7191 3453-211
E-Mail: info@sis.gmbh

For further information and documentation please refer to:
<https://www.soehnle-professional.com/en/site/documents>

This manual has been written for the following models:
► Veterinary scale 6958.01.002

Please read this manual completely before installation and operation.

1. Scope of delivery	26
2. Warnings	26
3. General information	27
3.1 Technical Data	27
3.2 Connection types.....	28
4. Installation	29
4.1 Unpacking and installation	29
5. Controls and functions	31
5.1 Display unit and key descriptions.....	31
5.2 Display character definitions	32
5.3 Terminal display.....	33
5.4 Meaning of symbols on faceplate.....	33
6. Operation	34
6.1 Normal weighing mode.....	34
6.2 Count weighing mode	36
6.3 Percent weighing mode.....	37
6.4 BMI working mode.....	39
6.5 HOLD function	40
6.6 Accumulation	41
6.7 Weight fine-tuning.....	42
7. Explanation of symbols & Troubleshooting.....	43
7.1 Explanation of symbols	43
7.2 Troubleshooting.....	44

1. Scope of delivery

- ▶ Terminal ✓
- ▶ Platform ✓
- ▶ User manual ✓

2. Warnings

- ▶ Verify that the input voltage range printed on the data label matches the local AC power to be used.
- ▶ Make sure that the power cord does not pose a potential obstacle or tripping hazard. (Power supply unit optionally available)
- ▶ Use only approved accessories and peripherals.
- ▶ Operate the equipment only under ambient conditions specified in these instructions.
- ▶ Disconnect the equipment from the power supply when cleaning.
- ▶ Do not operate the equipment in hazardous or unstable environments.
- ▶ Do not immerse the equipment in water or other liquids.
- ▶ Service should only be performed by authorized personnel.

3. General information

3.1 Technical data


Max. display range	-999,999 to 999,999
Max. display resolution	1:100,000 (primary unit) 1:125,000 (second unit)
Max. approved resolution	1:5,000
Division	0.0001, 0.0002, 0.0005 to 10, 20, 50
Calibration unit	kg or lb
Weighing units	g, lb, lb:oz, g, oz, pcs, %
Display	6-digit, 7-segment, 25 mm height LCD with blue backlight
Accuracy	≤ 0.01 %
Power supply	4x AA batteries or AC adapter: 9V DC / 600 mA with central positive
Work current	≤ 15 mA (backlight off, no load cells) ≤ 30 mA (backlight on, no load cells) ≤ 50 mA (backlight on, one load cell)
Functions	Normal weighing, counting, percent weighing, check weighing, BMI, Hold, Accumulation
Zero range	programmable zero range
Tare range	full capacity
Stabilization time	< 3 seconds
Keypad	[HOLD/SETUP], [PRINT/FUNC], [ACC/TOTAL], [UNIT/DATA], [TARE/PRESET], [ZERO/ON/OFF]
Operating temperature	-10° to 40°C
Humidity range	< 90% relative humidity, non-condensing
Load cell excitation voltage	5V DC
Load cell drive	-3mV/V to 3mV/V
Recommended actual scale interval (d)	> 1 uV / display division
Max. load cell excitation	60 mA
Signal connection	4 or 6 lead with sense leads
Real clock	built-in nonvolatile real time & date
Dimensions of scale	1.043,5 x 1.010 x 63 mm per platform resp. 2.087 x 1.010 x 63 mm complete scale
Dimensions of terminal	with bracket: 250 x 125 x 70 mm without bracket: 200 x 95 x 45 mm

3. General information


3.2 Connection types

7-pins socket, assigned as follows:

- ▶ Pin 1: Excitation + -
- ▶ Pin 2: Sense +
- ▶ Pin 3: Signal +
- ▶ Pin 4: Excitation -
- ▶ Pin 5: Sense -
- ▶ Pin 6: Signal -
- ▶ Pin 7: Shield

7 holes socket (Terminal)

7 pins socket (Base)


4. Installation

4.1 Unpacking and installation


Caution:

Do not pull on the load cell cable and connectors during installation, as this may damage the parts.

1. Remove the two platforms 1 & 2 with the same serial number from the boxes and align them as shown below. Remember to unplug the two connection cables and plugs on each platform..


2. Lift platform 1 as shown below and plug the connection cable of platform 2 into the connection socket on the terminal box of platform 1.


4. Installation

3. Mount the mounting brackets as shown in the picture below to align two platforms. Note: It is necessary that the two platforms are loosely connected.

4. Plug the connection cable of platform 1 into the connection socket on the display.


5. Controls and functions

5.1 Display unit and key descriptions

Key	Mode	Description	
HOLD/ SETUP	Weighing, counting or percent mode	< 3 seconds	Enters or exits HOLD mode
		> 3 seconds	Enters SETUP mode
	Input data mode	< 3 seconds	Returns to last sub-menu
		> 3 seconds	Input decimal point
	Menu selection mode		Returns to last sub-menu
PRINT/ FUNC	Weighing, counting or percent mode	< 3 seconds	Sends output data via the serial port
		> 3 seconds	Selects mode: Weighing, counting or percent
	Input data mode		Increases digit in flashing data entry position by one
	Menu selection mode		Returns to last item of current sub-menu
ACC/ TOTAL	Weighing, counting or percent mode	< 3 seconds	Adds accumulation values to memory, displays instances and totals
		> 3 seconds	Displays accumulation instances and totals
	Input data mode		Decreases the digit in the flashing data entry position by 1
	Menu selection mode		Goes to the next item of current sub-menu
UNIT/ DATA	Weighing mode	< 3 seconds	Changes weighing unit of measure
	Counting or percent mode	< 3 seconds	Enters the submenu to input piece weight for counting or to enter reference weight for percent-weighing
	Weighing, counting or percent mode	> 3 seconds	Enters the submenu to input the comparative data range for check-weighing
		> 3 seconds	Enters time or date setting mode
	Input data mode		Shifts the flashing data entry position from right to left
	Display ADC code		Select displaying code from no-filter, filter1, filter2
	Menu selection mode		Goes to next item of current sub-menu
TARE/ PRE- SET	Weighing, counting or percent mode	< 3 seconds	Tare the weight
		> 3 seconds	Enters pre-determined tare input mode
	Input data mode		Confirms the input data and forwards to next step
	Menu selection mode		Confirms the input data and forwards to next step
ZERO/ ON- OFF	Power Off		Power On
	Weighing, counting or percent mode	< 3 seconds	Zeros the platform weight
		> 3 seconds	Power off
	Input data mode		Ignores the modification
	Menu selection mode		Exits from current working mode


5. Controls and functions

5.2 Display character definitions

ASCII	LCD/LED Show	ASCII	LCD/LED Show	ASCII	LCD/LED Show
0		A		N	
1		B		O	
2		C		P	
3		D		Q	
4		E		R	
5		F		S	
6		G		T	
7		H		U	
8		I		V	
9		J		W	
		K		X	
		L		Y	
		M		Z	

5. Controls and functions

5.3 Terminal display

5.4 Meaning of symbols on faceplate

▶ Scale is zeroed, gross weight is 0, tare is 0

▶ Scale is stable

NET

▶ Display reading is net weight; tare is not 0

Total

▶ Display data is accumulated total times, weight, pieces, or percentage

Hold

▶ Scale is in dynamic weighing mode
 - Hold flashes - actual fluctuating weight is displayed
 - Hold does not flash - locked weight is displayed

Peak

▶ Peak - Scale is in dynamic weighing mode. Hold type is PEAK-HOLD

BMI

▶ Scale is working in BMI mode

lb

▶ Measure unit is lb or lb:oz

oz

▶ Measure unit is oz or lb:oz

kg

▶ Measure unit is kg

g

▶ Measure unit is g

%

▶ Measure unit is % (in percentage weighing mode)

Pcs

▶ Measure unit is pieces (in counting mode)

▶ Battery level

HI

▶ Data compare (check-weighing) is enabled. Current data (weight, pieces, or percent) is above the specified upper limit

OK

▶ Data compare is enabled. Current data is between the specified upper and lower limits

LO

▶ Data compare is enabled. Current data is below the specified lower limit

6. Operation

6.1 Normal weighing mode

Note:

During key operation, please note to use the second function of one key need pressing the key over 3 seconds. To input data or select menu, use arrow keys to process

▶ **Power on/off terminal:**

Power on: When terminal is off, press [ZERO/ON/OFF] key short to turn on.

Power off: When terminal is on, press [ZERO/ON/OFF] key long to turn off.

▶ **Zero:**

When the reading is stable and within the zero range, press [ZERO/ON/OFF] key to set new zero point. Please refer to the Zero and Tare limitations in the chart below.

▶ **Tare:**

When the gross weight is bigger than zero and the scale is stable, press [TARE/PRESET] key.

The terminal will show net weight of zero and the NET symbol will be lighted.

Please refer to the Zero and Tare limitations in the chart below.

▶ **Preset tare weight:**

Press [TARE/PRESET] key long, then „Pr.Tare“ will be shown and the „Tare“ symbol flashes. The scale is in Tare weight mode now. Use the arrow keys to input tare weight. The unit will be the same as before. There is no limitation to preset tare weight, but it should be bigger than zero.

After having input a tare weight, the NET symbol will be lighted.

Note:

This terminal can only save one tare weight. The new tare weight will automatically replace the old one. Please refer to the Zero and Tare limitations in the chart below.

▶ **Clear tare weight:**

Remove any weight from the platform, wait until the scale is stable, press [TARE/PRESET] key short. Please refer to the Zero and Tare limitations in the chart below.

▶ **Limitations to Zero and Tare operations under different conditions:**

Standard	Weight on platform	Data in Tare memory unit	Key function	
			Tare key	Zero key
Europe	≤ 0	No	No action	Zero
		Yes	Clear the tared weight	Zero and clear the tared weight
	> 0	No	Tare	Zero
		Yes		Zero and clear the tared weight

Note:

▶ ZERO is only active, when scale is stable and weight is in setting range.

▶ TARE is only active, when scale is stable.

▶ Clear TARE weight or ZERO scale, will make terminal to enter displaying GROSS mode.

▶ TARE weight, will maketerminal to enter displaying NET mode.

6. Operation

► **Change weight unit:**

Press [UNIT/DATA] key short, to select the units kg, lb, lb:oz, g or oz.

Note: Under some conditions, the units g and lb:oz are not available.

► **Check weight in weighing mode:**

- (1) To make weight compare function be available, *COmpARE-CORPAr* item should set to *YES*, and high and low limitation of weight should be set correctly according to following steps:
- (2) In weighing mode, press [UNIT/DATA] key for more than 3 seconds to input compare data of high and low. When *Unit L_H9* is shown, use [PRINT/FUNC] or [ACC/TOTAL] key to select the unit of weight. Use [TARE/PRESET] key to confirm and go to next. Press [ZERO/ON/OFF] key to exit and go back to weighing mode.
- (3) After *Hi 9H* being shown, *000000* will be displayed, use the arrow keys to input high weight number and press [TARE/PRESET] key to confirm. The symbol „Hi“ will be shown in this step. Press [ZERO/ON/OFF] key to exit and go back to weighing mode.
- (4) After *Lo 1'* being shown, *000000* will be displayed, use keys to input low weight number and press [TARE/PRESET] key to confirm. The symbol „Lo“ will be shown in this step. Press [ZERO/ON/OFF] key to exit and go back to weighing mode.
Note: If High number is 0 or is equal or less than low number, the comparison will be disabled, and the input data has no limitation.
- (5) After a reasonable limitation is set and compare is be active, one of the symbols HI, OK, LO will be lighted, and the beeper will sound according to its setting in *USER-beEP*.

► **Accumulation:**

Press [ACC/TOTAL] key to add displayed number to accumulation memories, and accumulation times will also add up 1. And then to display accumulation result.

► **Data output:**

When reading is stable, press the [PRINT/FUNC] key.

► **Change working mode:**

Press [PRINT/FUNC] keylong, then use the arrow keys to choose and confirm to enter into weighing mode, counting mode, percentage working mode or BMI working mode.

► **Enter to or exit from HOLD mode:**

Press the [HOLD/SETUP] key.

► **Enter Setup mode:**

- (1) If you need to set configuration parameters, set user parameters, calibrate the scale, set current date or time, test some hardware. It is necessary to press the [HOLD/SETUP] key long, to enter setup mode
- (2) After entering Setup mode, the main menu item *COmpARE* will be shown first.
- (3) In Setup mode, use the arrow keys to select the wanted submenu and the wanted menu item, select wanted choice, set wanted number, confirm and save data or exit this mode.

6. Operation

6.2 Count weighing mode

- ▶ In this mode, scale will weigh goods weight on scale, calculate and display its counts after the piece weight of goods is obtained.
- ▶ To make counting function be available, `COUING-FUNC-COUNT` item should be set to `YES` in `COUING` menu.
- ▶ **Enter counting working mode:** In normal weighing or percent-weighing mode, press [PRINT/FUNC] key long. `WEIGH` or `PERCENT` or `UNIT` will be shown, use [ACC/TOTAL] or [PRINT/FUNC] key to select `COUNT`. Then press [TARE/PRESET] key to confirm change to parts counting mode. The weight of the last determined piece will be used until the weight of a new piece is taken over.

In counting working mode, the functions ZERO, ON/OFF, TARE, PRESET TARE, ACC, PRINT, HOLD, SETUP are available.

- ▶ **There are two ways to obtain the piece weight:**
 - (1) Direct weight of input piece
 - (2) Weight of samples with known quantity
- ▶ **Input piece weight from keypad:** In counting mode, press the [UNIT/DATA] key, when `INP_WT` is shown, press [TARE/PRESET] key to enter input piece weight mode.
 - (1) When `UNIT_H` is shown, use [PRINT/FUNC] or [ACC/TOTAL] key to select the unit of piece weight, use [TARE/PRESET] key to confirm and go to next. Press [ZERO/ON/OFF] key to exit getting piece weight mode and back to counting mode.
 - (2) When last stored piece weight is shown, use the arrow keys to input the new piece weight. Press [HOLD/SETUP] key for more than 3 seconds to input a decimal point. Press [TARE/PRESET] key to confirm and save it. Then return to counting mode. If the input piece weight is less than 0.5d, the terminal will display `PLT_Er` and returns to counting mode.
- ▶ **Obtain piece weight by weighing samples weight which quantity is known:** In counting mode, press [UNIT/DATA] key, When `INP_PCS` is shown, use [PRINT/FUNC] or [ACC/TOTAL] key to select `SPL_PCS`, press [TARE/PRESET] key to weigh samples weight (which quantity is known), then calculate the piece weight. Press [ZERO/ON/OFF] key to exit getting piece weight mode and return to counting mode.
 - (1) When `SPL_Lb` is shown, move away any sample on scale and press [TARE/PRESET] key to confirm, before scale is stable, `SPL_Lb` will be flashed. After it is stable, it will go to the next step. Press [ZERO/ON/OFF] key to exit getting piece weight mode and back to counting mode.
 - (2) When `SPL_H` is shown, put samples (its quantity is known) onto the scale, Press [TARE/PRESET] key to confirm reading weight. Before scale is stable, `SPL_H` will be flashed. After it is stable, it will go to next step. Press [ZERO/ON/OFF] key to exit getting piece weight mode and back to counting mode.
- ▶ After `INP_PCS` being shown, `000001` will be displayed, use the arrow keys to input the quantity of samples and Press [TARE/PRESET] key to confirm. If the calculated piece weight is less than 0.5 d, the terminal will display `PLT_Er` and go back to counting mode, otherwise, after the reasonable piece weight being got, the scale will go back to counting mode. The gotten piece weight can be saved after the power off and can be used next time.

6. Operation

► **Check Counts (counts compare) in Counting mode:**

- (1) To make counts compare function be available, `COU/FI 9-FUNC-COMP` item should set to `YES`, and high and low limitation of pieces should be set correctly according to following steps:
- (2) In counting working mode, press [UNIT/DATA] key for more than 3 seconds to input compare data of high and low.
- (3) After `HI 9H` being shown, last upper limitation will be displayed, use the arrow keys to input high quantity number and press [TARE/PRESET] key to confirm. Symbol `Hi` will be shown in this step. Press the [ZERO/ON/OFF] key to exit getting piece weight mode and back to counting mode.
- (4) After `Lo 9L` being shown, last lower limitation will be displayed, use the arrow keys to input low quantity number and press [TARE/PRESET] key to confirm. The symbol of `Lo` will be shown in this step. Press [ZERO/ON/OFF] key to exit getting piece weight mode and back to counting mode.
Note: If the high number is 0 or is equal or less than low number, the comparison will be disabled.
- (5) After a reasonable limitation is set and compare is be active, one of annunciators `HI`, `OK`, `LO` will be lighted, and the beeper will sound according to its setting in `USER-BEEP`.

6.3 Percent weighing mode

- In this mode, scale will weigh goods weight on it, calculate and display its percentage after the unit-percentage-weight of goods is obtained.
Note: If 100% display format is set to 100%, 100.0% or 100.00% in `COU/FI 9-FUNC-PERCEN` menu item, then, the unit-percentage-weight is the weight of 1%, 0.1% or 0.01%.
- To make percent weighing function be available, `COU/FI 9-FUNC-PERCEN` menu item shouldn't be set to `no`.
- To enter percent weighing mode, in normal weighing or counting mode, press the [PRINT/FUNC] key long, `WEI 9H` or `COU 9L` or `bit` will be shown, use [ACC/TOTAL] or [PRINT/FUNC] key to select `PERCEN`, then press [TARE/PRESET] to confirm go to percent weighing mode.
Before new unit-percentage-weight is available, the last unit-percentage-weight will be used.
The last unit-percentage weight will be used until a new unit-percentage weight is taken over.
- In percent weighing mode, the functions ZERO, ON/OFF, TARE, PRESET TARE, ACC, PRINT, HOLD, SETUP are available.
- To obtain the unit-percentage-weight, there're two ways:
 - (1) Input weight and its percentage, then scale calculates the unit-percentage-weight.
 - (2) Weigh samples weight which percentage is known

6. Operation

- ▶ Input weight and its percentage from keypad, and calculate unit-percentage-weight: in percent weighing mode, press [UNIT/DATA] key, When $1nP_Pct$ is shown, press [TARE/PRESET] key to enter this mode:
 - (1) Before input weight, use [PRINT/FUNC] or [ACC/TOTAL] key to select the percentage from 1%, 2%, 5%, 10%, 20%, 50% and 100%, this percentage is corresponding to the weight you will input in the following steps
 - (2) When $Unit_H9$ is shown, use [PRINT/FUNC] or [ACC/TOTAL] key to select the unit of input weight, use [TARE/PRESET] key to confirm and go to next. Press [ZERO/ON/OFF] key to exit and return to percent weighing mode.
 - (3) When last stored unit-percentage-weight data is shown, use keys to input new unit-percentage-weight, press [HOLD/SETUP] key more than 3s to input decimal point. Press [TARE/PRESET] key to confirm and save it, then go back to percent weighing mode. If the calculated unit-percentage-weight is less than 0.5d, the terminal will display Pct_Er and return to percent weighing mode.

- ▶ Obtain unit-percentage-weight by weighing samples weight which percentage is known: in percent weighing mode, press [UNIT/DATA] key, When $1nP_Pct$ is shown, use [PRINT/FUNC] or [ACC/TOTAL] key to select $5PL_Pct$, press [TARE/PRESET] key to weigh samples (which percentage is known) weight, calculate piece weight.
Press [ZERO/ON/OFF] key to exit and return to percent weighing mode.
 - (1) When $5PL_Lo$ is shown, remove any sample on scale and press [TARE/PRESET] key to confirm, before scale is stable, $5PL_Lo$ will be flashed. After it is stable, it will go to next step.
Press [ZERO/ON/OFF] key to exit and return to percent weighing mode.
 - (2) When $5PL_Hi$ is shown, put samples (with known percentage) onto the scale.
Press [TARE/PRESET] key to confirm reading weight. Before scale is stable, $5PL_Hi$ will be flashed. After it is stable, it will go to next step. Press [ZERO/ON/OFF] key to exit and back to percent weighing mode.
 - (3) After $1nP_Pct$ being shown, 000000 (position of decimal point is determined by $COVF9-FUNC-PErCEr$ setting) will be displayed. Use the arrow keys to input the percentage of samples.
Press [TARE/PRESET] key to confirm. If the calculated unit-percentage-weight is less than 0.5 d, the terminal will display Pct_Er . Returns to percent weighing mode, otherwise, after the reasonable unit-percentage-weight being got, the scale will go back to percent weighing mode. The gotten unit-percentage-weight can be saved after the power off and can be used next time.

- ▶ Check Percent (percentage compare) in Percent weighing mode:
 - (1) To make percentage compare function be available, the $COVF9-FUNC-COvPAR$ menu item should be set to YES and high and low limitation of percentage should be set according to the following steps.
 - (2) In percent weighing mode, press [UNIT/DATA] key for more than 3 seconds to input compare data of high and low.
 - (3) After $Hi 9H$ being shown, 000000 will be displayed. Use the arrow keys to input high percentage number and press [TARE/PRESET] key to confirm. The symbol of Hi will be shown in this step.
Press the [ZERO/ON/OFF] key to exit and return to percent weighing mode.
 - (4) After $Lo 1'$ being shown, 000000 will be displayed. Use the arrow keys to input low percentage number and press [TARE/PRESET] key to confirm. The symbol of Lo will be shown in this step.
Press [ZERO/ON/OFF] key to exit and back to counting mode.
Note: If High number is 0 or is equal or less than low number, the comparison will be disabled.
 - (5) After a reasonable limitation is set and compare is active, one of the symbols HI, OK, LO will be lighted, and the beeper will sound according to its setting in $USER-bEEP$.

6. Operation

6.4 BMI working mode

- ▶ To make BMI working Mode be available, **LOAD-FUNC-BMI** menu item should be set to **YES** and factory setting should be also enable this function.
- ▶ To enter BMI working mode:
 - ▶ When **LOAD-FUNC-ACCUTIME=YES**: If In normal weighing mode, percent weighing mode, or counting mode, press [PRINT/FUNC] key long, then **WEIGHT** or **LOAD** or **PERCENT** will be shown, use [ACC/TOTAL] or [PRINT/FUNC] key to select BMI, then press [TARE/PRESET] to confirm go to BMI mode.
 - ▶ When **LOAD-FUNC-ACCUTIME=NO**: press [ACC/TOTAL] (BMI) key, to go to BMI mode.
 - ▶ After scale went to this mode, "**HT.xxx**" (means: last input height is xxx cm) will be displayed. Wait for input height:
 - (1) To change height unit to cm or inch by pressing [UNIT/DATA] key
 - (2) To change height number by using [PRINT/FUNC], [ACC/TOTAL] (BMI) keys
 - (3) Press and hold [PRINT/FUNC] or [ACC/TOTAL] (BMI) key will increase or decrease number
 - (4) Press [TARE/PRESET] key to confirm the input.
Press [ZERO/ON/OFF] key to exit input data mode and return to BMI working mode.
The range of height is 50-250cm and default is 170cm.
 - ▶ In this mode, when BMI number is shown (BMI symbol is also on), or weight number is shown (BMI and kg or lb annunciators are on), Press [ACC/TOTAL] key to select weight or BMI number to be displayed, when weight is displayed, the weight unit can be selected by pressing [UNIT/DATA] key, and BMI and weight unit will be displayed at same time.
 - ▶ In this mode, when current net weight is less than **LOAD-RNG**, the indicator will go to display weight number if **LOAD-FUNC-ACCUTIME=NO**; or the indicator will back to original working mode if **LOAD-FUNC-ACCUTIME= YES**.

6. Operation

6.5 HOLD function

- ▶ HOLD function can be used to freeze display number. In this mode, the scale can catch a dynamic number, hold a stable number, or average a unstable number, then HOLD (freeze) this number temporary for user to watch or record. This function can be used in normal weighing mode, counting mode and percent weighing mode. The speed of A/D converter is increased to 80Hz (if *CONF19-Ad_H_Spd* is set to *YES*) from original 10Hz for some dynamic weighing applications. With the hold function, it is possible to weigh restless weighing samples such as live animals, moving objects. The indicator provides special mode settings to accommodate sample's movements.

- ▶ To make HOLD function be active, the *CONF19-FUNC-HOLD* menu item must be set to *YES*; menu items of *USER-HOLD-HLD_00d /-Rd9_00d /-HLD_00d /-dYn_rn9 /-Stb_00d*, *USER-0tHEr_nLd_rn9* need to be set at a reasonable value. To speedup sampling of weight, set *CONF19-Ad_H_Spd* menu item to *YES*. Please note: lower speed will get more stable data, but some instant weight information may be lost. To enter HOLD working mode, press [HOLD/SETUP] key when scale works in normal weighing mode, counting mode or percent weighing mode.

- ▶ There are several HOLD mode to freeze display data:
 - (1) Positive Peak Number HOLD mode
 - (2) Negative Peak Number HOLD mode
 - (3) Toggle HOLD mode
 - (4) Average HOLD mode
 - (5) Auto HOLD mode

- ▶ **Positive Peak HOLD:**

When *USER-HOLD-HLD_00d* is set to *PS_PEAH*, the hold mode is positive peak hold mode. When scale first enters this working mode, it will display the largest positive number that is from the time of zero-point set. After entering this working mode, scale will always catches and refresh positive larger number and display it. To exit HOLD mode, press [HOLD/SETUP] key again.

- ▶ **Negative Peak HOLD:**

When *USER-HOLD-HLD_00d* is set to *ng_PEAH*, the hold mode is negative peak hold mode. When scale first enters this working mode, it will display the largest negative number that is from the time of zero-point set. After entering this working mode, scale will always catches negative larger number and display it. To exit HOLD mode, press [HOLD/SETUP] key again.

- ▶ **Toggle HOLD:**

When *USER-HOLD-HLD_00d* is set to *0099LE*, the hold mode is toggle hold mode, a manual Hold function. After entering this working mode, the scale will freeze and display the number if scale is stable. Only the weight that is over *USER-0tHEr_nLd_rn9* (zero ,dead' band) can be held. To exit HOLD mode, press [HOLD/SETUP] key again. If the time of scale being unstable is more than *USER-HOLD-5tb_00d*, *5tb_Er* will be shown, press [TARE/PRESET] key to start averaging again or press [HOLD/SETUP] key to exit.

6. Operation

► **Average HOLD:**

When *USER-HOLD-HLD-TOd* is set to *AUERAG*, the hold mode is average hold mode. After entering this working mode, scale will freeze and display number if scale is stable. If scale is not stable, but the variation is less than *USER-HOLD-dyn-rn9*, scale will average data in *USER-HOLD-Au9-ti7*, then freeze and display the number. Only the weight that is over *USER-DEHER-nLd-rn9* can be frozen. Scale will exit HOLD mode according to the setting of *USER-HOLD-HLD-ti7*. If the time of scale variation being over *USER-DEHER-nLd-rn9* is more than *USER-HOLD-5tb-ti7*, *5tb-Er* will be shown, press arrow keys to start averaging again, or press HOLD/SETUP key to exit.

► **Auto HOLD:**

When *USER-HOLD-HLD-TOd* is set to *AUTd*, the hold mode is a auto hold mode. Different subjects can be weighed one after another without pressing any buttons. After entering this working mode, scale will freeze and display number if scale is stable. If scale is not stable, but the variation is less than *USER-HOLD-dyn-rn9*, scale will average data in *USER-HOLD-Au9-ti7*, then freeze and display the number. Only the weight that is over *USER-DEHER-nLd-rn9* can be frozen. If the held weight is moved away, and a new load put on the scale, scale will automatically hold new number of load. Scale will exit HOLD mode according to the setting of *USER-HOLD-HLD-ti7*. If the time of scale variation being over *USER-DEHER-nLd-rn9* is more than *USER-HOLD-5tb-ti7*, *5tb-Er* will be shown, press [TARE/PRESET] to start averaging again, or press [HOLD/SETUP] key to exit.

► In positive or negative Peak HOLD mode, the PEAK and HOLD symbol will be lighted. In another HOLD mode, the HOLD symbol will be lighted. When the HOLD symbol flashes, the displayed number is active. When the HOLD symbol becomes steady, the displayed number is frozen.

6.6 Accumulation

- Accumulation function can be used in normal weighing mode, counting mode and percent weighing mode, When this function is enabled, you can accumulate current net weight, piece, and percentage. Note, only the load on scale is larger than *USER-DEHER-nLd-rn9*, the displayed positive number can be added up. The accumulation times and total can be displayed or printed.
- To make data accumulation function be available, *COntrol-FUNC-ACCUM* menu item should set to *normal* or *AUTd*.
- When *COntrol-FUNC-ACCUM* is set to *normal*, the stable and positive displayed net weight (must be larger than *USER-DEHER-nLd-rn9*), piece or percentage can be accumulated by long pressing [ACC/TOTAL] key, and the terminal will display accumulation times first. Then display total of number. To avoid repeating accumulation for same load, one load only can be accumulated once. So, before a new load put onto the scale, the original load should be removed and let load on scale be smaller than *USER-DEHER-nLd-rn9*.
- When *COntrol-FUNC-ACCUM* is set to *AUTd*, the stable and positive displayed net weight (must be larger than *USER-DEHER-nLd-rn9*), piece or percentage can be accumulated automatically, and indicator will display accumulation times first, and then display total of number. To avoid repeating accumulation for same load, one load only can be accumulated once. So, before a new load put onto the scale, the original load should be removed and let load on scale be smaller than *USER-DEHER-nLd-rn9*.
- To view total, when display number is zero, press [ACC/TOTAL] key long. It will alternatively display the accumulation times and the accumulated sum total thus far (weight or quantity), until the [ACC/TOTAL] key being pressed again. Note: When HOLD function is enable, and scale is working in PEAK HOLD mode (*COntrol-HOLD-ZYES*, *USER-HOLD-HLD-TOd=PS-PEAK-rn9-PEAK*), accumulation function will be automatically disabled!

6. Operation

6.7 Weight fine-tuning

- ▶ With this function, the user can adjust the displayed weight a little. There is no need for a standard weight. This adjustment can only be done in the weighing mode.

Note:

- (1) The scale must have been calibrated before this adjustment
 - (2) The range of adjustment is "(current displayed weight) x (0.9-1.1)".
That means the range is about $\pm 10\%$
 - (3) The "CONF:G-FUNC-LE_ADJ=YES" must be set.
 - (4) Only an authorized manufacturer's representative or certified verification personnel are authorized to make these changes. Changing this value alters the calibration values!
- ▶ To enter this mode, turn on the terminal. After the terminal is displaying 0 weight, put a load onto the scale (e.g. the correct weight is 10.0 kg), then the terminal will display the load's weight, saying "10.5 kg". Press [TARE/PRESET] and [ZERO/ON/OFF] keys at the same time until the first digit flashes. This means indicator has entered into "weight fine-tune" mode.
 - ▶ Use the arrow keys to input the correct weight (10.0). After having confirmed by [TARE/PRESET], the active correct weight will be displayed and no any digit will be flashed. After this, the displayed weight will be adjusted by this ratio (10.0/10.5) and this ratio will keep active until the next modification on it.
 - ▶ Procedure for the removal of this ratio's effect:
Move away weight on scale, press [ZERO/ON/OFF] to let 0 be displayed. Put a load onto the scale. A number will be displayed, suppose it's 10.0 kg. Press [TARE/PRESET] and [ZERO/ON/OFF] at same time until first digit flashes. This means the terminal has entered into "weight fine-tune" mode. Press [HOLD/SETUP] key. The displayed weight will be restored to 10.5 kg. Then press [TARE/PRESET] to confirm and exit to the normal weighing mode.

7. Explanation of symbols & Troubleshooting

7.1 Explanation of symbols

<i>0----</i>	Zero is over the setting range
<i>0....</i>	Zero point is below the setting range
<i>Ad-----</i>	Signal to ADC is over max. range
<i>Ad....</i>	Signal to ADC is below min. range
<i>----</i>	Weight is over upper limitation, or display data is over limitation
<i>----</i>	Weight is below lower limitation
<i>EEP.E1</i>	CONFIG or CAL parameters are not correctly set
<i>EEP.E2</i>	USER parameter is not correctly set
<i>Lo.bAt</i>	Battery voltage is lower than setting.
<i>CAP - - -</i>	Next displaying content is Capacity
<i>CAP.Er</i>	Parameters about Capacity is not correct
<i>CAL.PH</i>	Calibration on point(x)
<i>CAL.OFF</i>	Calibration Seal Switch is on OFF position
<i>CAL.ON</i>	Calibration Seal Switch is on ON position
<i>CAL.Er</i>	Calibration error, maybe input data or loaded weight is incorrect, unstable, un-linear
<i>CAL.End</i>	Calibration is end
<i>oFF</i>	Terminal is turned off
<i>Stb.Er</i>	Unstable time is lager than setting of USER-HOLD-STB.TIM
<i>RcC.HHH</i>	Accumulation times is xxx
<i>Pr.tArE</i>	Preset TARE weight
<i>COmP</i>	Go to COMPARE data mode
<i>HI9H</i>	Input HIGH limitation data of comparison
<i>LO!</i>	Input LOW limitation data of Comparison
<i>SPL.Lo</i>	Sample load weight of low point.
<i>SPL.Hi</i>	Sample load weight of high point.
<i>SPL.Pwt</i>	Sample goods weight to calculate piece weight
<i>InP.PCS</i>	Input pieces number of weighted goods
<i>Unit.kg</i>	Unit kg is selected
<i>Unit.lb</i>	Unit lb is selected
<i>Pwt.Er</i>	Piece weight error, too small (<0.5d)
<i>SPL.PCt</i>	Sample goods weight to calculate
<i>InP.PCt</i>	Input percentage of weighted goods
<i>PCt.Er</i>	Unit-Percentage -Weight is too small (<0.5d).
<i>RCU.Er</i>	Internal calculation overflows

7. Explanation of symbols and Troubleshooting

7.2 Troubleshooting

Problem	Probable cause	Remedy
<i>Ad---</i>	Load cell wires to indicator are incorrectly connected, or shorted, or opened; or ADC, load cell are damaged	Make sure wires are ok and correctly connected. Replace load cell or ADC chip. Service is required.
<i>Ad---</i>		
<i>0----</i>	Weight reading exceeds Power On Zero limit.	Make sure scale platform is empty. Perform zero calibration.
<i>0-----</i>	Weight reading below Power On Zero limit.	Install platform on scale. Perform zero calibration.
<i>-----</i>	Weight reading exceeds Overload limit, or the weight value cannot be displayed in the current unit of measure because it exceeds 6 digits.	Reduce load on scale until weight value can be displayed. Use a more appropriate unit of measure. Re-set some parameters of CONFIG or UAER .
<i>-----</i>	Weight reading below under load limit.	Install platform on scale. Perform zero calibration.
<i>EEP.E1</i>	CONFIG or CAL parameters are not correctly set	Reset items in CONFIG mode, do calibration
<i>EEP.E2</i>	USER parameter is not correctly set	Reset items in USER mode
<i>CAP.Er</i>	Capacity parameters are not correct	Set PRIM.N/PRIM.d/SECND.n to correct number, make sure capacity not more than 6 digit
<i>CAL.Er</i>	Calibration error, maybe input data or loaded weight is too small, too big, unstable, un-linear	Input correct data, load correct weight onto platform. Service is required.
<i>Pwt.Er</i>	Piece weight is error, it's too small (<0.5d). The weight on the platform is too small to define a valid reference weight.	Use a greater weight for the sample.
<i>Pct.Er</i>	Unit-Percentage -Weight error, too small (the weight of 1%, 0.1%, or 0.01% determined by CONFIG-FUNC-PERCNT is less than 0.5d)	Use a greater weight for the sample.
<i>CALC.Er</i>	Internal calculation overflow	Adjust the value of the PWT or PCT
<i>Stb.Er</i>	USER-HOLD-STB.TIM is too short. / USER-HOLD-HLD.RNG is too small. / Other failure	Set USER-HOLD-STB.TIM longer, or set USER-HOLD-HLD.RNG bigger. / Service is required
Does not turn on	Power cord not plugged in or properly connected. Power outlet not supplying electricity. Battery discharged. / Other failure.	Check power cord connections. / Make sure power cord is plugged into the power outlet. / Check power source. / Replace batteries. / Service is required.
Cannot zero display or will not zero when turned on.	Load on scale exceeds allowable limits. Load on scale is not stable. Load cell damage.	Remove load on scale. Wait for load to become stable. Service is required.
Cannot display weight in req. weighing unit.	Unit not set to enable, or $\geq 5\text{oz}$, when unit is lb:oz	Enable unit in CONFIG-UNITS
Battery symbol is empty or Lo.BAt is shown	Batteries are discharged.	Charge batteries

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All rights reserved.
© Soehnle Industrial Solutions GmbH, Publication, duplication, or any form of commercialization of such material beyond the scope of the copyright law shall require the prior written consent of Soehnle Industrial Solutions GmbH. Subject to technical modifications.